

SPARKY



The Sparky is designed to deliver full rehab capability in a compact, powerful paediatric package. Engineered to operate in tight child-sized spaces and to keep up with today's active child, the Sparky is an exceptional fit from the playground to the classroom and beyond.



CRASH
TESTED

Jazzy[®]
POWER CHAIRS



Distributed by:

MMS
MEDICAL



@PRIDEMOBILITYUK

CORK Tel: 021 461 8000 | Fax: 021 461 8099 | DUBLIN Tel: 01 450 0582 Fax: 01 456 5889

Email: info@mmsmedical.ie | Web: www.mmsmedical.ie



Coming Soon...

Sparky Plus



Features:

- 14" and 16" contoured back
- Stealth Thoracics
- iFit comfort harness
- iFit pelvic support

Jazzy Sparky | Colour Options & Tech Spec

Seating and Power Options

Ion Seat Power / Manual Options:

The Sparky comes a standard with Ion Power Tilt / Attendant Recline

Ion Seat Sizes: 10" x 10", 12" x 12" and 14" x 14"

(COMING SOON)

ADI Sizes: 14" and 16" Contoured back

Standard Features



Power Tilt



Attendant Recline

Colour Options



White



Yellow

Technical Specifications

Drive Wheels (cm) 30.48 cm pneumatic	Caster Wheels (cm) Front and Rear: 20.32 cm	Max Speed (mph) 4 Mph	Ground Clearance (cm) Battery Box: 9.5 cm	Turning Radius (cm) 74.3 cm (29.25")	Overall Length (cm) 95.9 cm (37.75")	Overall Width (cm) 61 cm (24")	Seating Options (inches) Ion Seating	Batteries x 2 40 AH Deep Cycle	Available Electronics VR2	Weight Capacity (stone) 8 Stone 13 lbs. (56.7 k.g)
Seat to ground height (") 16" - 16.5"	Base Weight (kg) 59.2 kg	Battery Weight Each (kg) 10.5 kg	Warranty <ol style="list-style-type: none">1. Lifetime limited warranty on frame.2. 1-year limited warranty on electronics.3. 18-month warranty on drive control system.4. 2-year limited warranty on drivetrain.				Notes <ol style="list-style-type: none">1. Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.3. Battery weight may vary based on manufacturer.4. Not available with iLevel Seat lift			

Get Social with Quantum



FOLLOW THE LIVES OF OUR QUANTUM STARS, AND GET THEIR PERSPECTIVE OF QUANTUM

www.quantumrehab.co.uk/quantumstars



@QUANTUMREHABUK

Ensure you have the most up to date information about Quantum Rehab by subscribing to our social media pages

Distributed by:



MMS Medical
51 East Gate Drive,
Little Island,
Cork.

MMS Medical
Unit 81
Western Parkway Business Park,
Lower Ballymount Road,
Dublin 12

Further information can be obtained through our web site, or prescription forms. All information is subject to change without notification. Please consult Quantum with any queries you may have.