The ultimate mobility tool

Batec handbikes are a comprehensive solution to your mobility problems:

- ✓ Go from reduced to **augmented mobility**; from needing help to being **independent** and **proactive**.
- ✓ Get a **2 in 1**, without giving up your manual wheelchair.
- ✓ Overcome obstacles and drive over bumpy terrains.
- ✓ Attach our handbikes in just seconds and with no need for transfers.
- Also available for quadriplegics and hemiplegics.



Request a quote or a demonstration without any commitment!



You CAN get used to being in a wheelchair

But not to giving up on

www.batec-mobility.com

BATEC MOBILITY S.L.

Illa de Buda, 2 08192 Sant Quirze del Vallès (Barcelona) - ESPAÑA 93 544 2003 <u>info@b</u>atec-mobility.com



BATEC ELECTRIC 2, the perfect balance

The **BATEC ELECTRIC** is the latest in our saga of handbikes. It features the latest **Batec Mobility** technology and an **18" Wheel** designed exclusively for this model, which means you can use it both in **urban** and **off-road settings**.





All-in-one LCD screen with gear selector.



18" Wheel, ultra-grip tyre and aluminium double-wall rim.



250 RPM 900W high torque brushless motor with reverse gear. Maximum speed



700-lumen front light with position lights / high bearn options. Dual LED tall light son



the stand frame.



Horn and USB port.



180 mm disc, two mechanical brakes,

two aluminium levers with parking brake.

Removable right dropout for easy inner tube or tyre change.



Removable QR weights.



21,5 km* autonomy with 36 V - 280 Wh battery. 44 km* autonomy with 36 V - 576 Wh battery. 4A charger.



Matte black aluminium handlebar with folding system for convenience.



G&B handlebar option without levers for quadriplegics.



Handlebar option with right or left Monolever brakes for hemiplegics.



^{*} Autonomy calculated under the following conditions: 75 - kg user, ultra-light wheelchair, flat terrain, no wind and steady speed of 15 km/h. May vary +- 20% depending on the battery. Number/version: MM-008/7.0. | October 2023