

Access Padded

The Oxford® Access Padded sling has been designed specifically to facilitate the toileting process by allowing excellent access to clothing. It is an easy to fit sling, fully padded for additional comfort. It is not a general purpose sling and will only suit around 25% of clients.

- ✓ Available with or without head support
- ✓ Adjustable body belt for added support and security
- ✓ Non-slip back section ensures the sling does not 'ride up'
- ✓ Colour coded binding indicates sling size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Can be customised to fulfil patient specific needs



WITH HEAD SUPPORT




Chair to Chair



Bed to Chair/Chair to Bed
(From seated position only)



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size	Part No. (Padded)	Part No. (Padded HS)	SWL
XS	SL1640	SL1625	227 kg
S	SL1641	SL1626	227 kg
M	SL1642	SL1627	227 kg
L	SL1643	SL1628	227 kg
XL	SL1644	SL1629	227 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.

