

POWER RECLINERS

C1 - Petite




DUAL MOTOR LIFT CHAIR



- Dual motor Layflat liftchair
- Petite lift chair
- Weight capacity of 25 Stone (158.7 kg)

Fabric



Burgundy Vinyl

Lift Chair Model	Back Width	Seat to Floor	Seat to Depth	Seat Width	Top of Back to Seat	Weight Capacity
C1	58.42 cm (23")	44.45 cm (17.5")	44.45 cm (17.5")	46 cm (18")	58.42 cm (23")	25 Stone (158.7 kg)

MMS Medical, 51 East gate Drive, Little Island, Cork. Tel : 021 4618000 Email: info@mmsmedical.ie

MMS Medical, Unit 81 Western Parkway Business Park, Lower Ballymount Road, Dublin 12. Tel : 01 4500582

www.mmsmedical.ie