

zippie®

TS™

Paediatric Folding Tilt and Positioning chair



Featured With New Glow in the dark paint

**MMS**  
MEDICAL

improving people's lives

The Zippie TS is the original premier paediatric tilt-in-space chair that is available in either a folding or rigid style frame. The Zippie TS has built-in frame growth with the widest variety of options and is able to accept many seating and positioning solutions.



Tilt-in-Space up to 45 degrees. Available in both rigid and folding frames. Also available with Foot Release Tilt Actuator, which eliminates exposed cables.

Ability to grow the chair frame as the child's needs change. The Zippie TS™ provides growth in both seat depth and seat width.

Optional Folding Frame is a lightweight welded frame construction which allows the frame to compress into a small package that is ideal for transportability and storage in tight places.



ZIPPIE TS SPECIFICATIONS	Folding w/Tilt	Rigid w/Tilt
Overall weight:*	13kg / 29 lbs. w/o footrest	13kg / 29 lbs. w/o footrest
Weight capacity:	75kg / 12 st.	75kg / 12 st.
Overall product width:	46 - 62cm (18.25" -24.25")	46 - 62cm (18.25" -24.25")
Transit option:**	Yes	Yes
Front seat height:	40 - 50cm (15.5" - 20")	40 - 50cm (15.5" - 20")
Rear seat height:	35.5 - 47cm (14" - 18.5")	35.5 - 47cm (14" - 18.5")
Hanger options:	60° S-A, 70°, 80°, 90° (S-A, Front Mount & HD), ELR, ALR	60° S-A, 70°, 80°, 90° (S-A, Front Mount & HD), ELR, ALR
True Fit Growth System:	Yes	Yes
Warranty:	Lifetime on frame & cross brace, 1 year on parts	Lifetime on frame & cross brace, 1 year on parts

\* Actual measurements may vary according to the wheelchair's configuration and accessories

\*\* Transit approved to ISO7176 Part 19 when the transit option is selected.

For comprehensive information on the features and functions of the Zippie TS please consult the Zippie TS Owners Manual. The Owners Manual can be downloaded from the Sunrise Medical website [www.sunrisemedical.com](http://www.sunrisemedical.com).

The Zippie TS Owners Manual can also be made available upon request or by contacting your nearest Sunrise Medical product provider.