

PRODUCT DATASHEET

Disposable Bedpan and Support, 2 L

Brand	: CURAS®
Category	: Nursing
Range	: Bedside Toileting Supports
Product REF	: 07000
Product Description	: Disposable Bedpan & Support, 2L, Blue, Large, Medium Duty (50-75kg), Single Patient Use (for REF 06305/06, 06403, 08000/09 & 09700), CleanPak™



PRODUCT DESCRIPTION

Disposable plastic bedpan with capacity of 2 L and can be used as a support.

INTENDED USE

It is intended to use to hold and transport human waste for safe disposal. A disposable bedside toileting support product designed specifically for single-patient-use on a short-term hospitalization. It is also designed to fit pulp bedpan liner (08000, 08009). It can be used together with Hygiene Control Bag™ (06305 and 06306) for better hygiene.

FEATURES & BENEFITS

- Designed specifically for single-use that reduce risks of cross-contamination and hospital-acquired infections.
- Economical & environmentally friendly product made of 7% recycled material.
- Durable, heavy duty plastic gives extra strong support for individuals weighing up to 115kg.
- Fits well on all standard commode chairs.
- Smooth surface for easy cleaning and disinfection.
- Soft surface for less traumatic patient contact.
- No sharp edges to prevent any kind of injuries.
- Stackable to allow it to be stored with great efficiency.
- CleanPak™ inner bag protects product from dust and contamination before reaching its point of use.

PRODUCT SPECIFICATIONS

REF	Dimensions (L x W x H)	Average Weight/Pc	Packing
07000	37.0 x 30.0 x 10.5 cm	155 g	CleanPak™

Composition

- The bedpan is made of polyethylene (PE).

Country of Origin

- China

Sterilization

- Non-sterile

Shelf Life

- 10 years

Selected Standards & Guidelines

- ISO 13485:2016

Recommended Storage Condition

- Store the product at room temperature and avoid exposure to extreme temperatures.
- Keep the product at a cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- If regular waste, dispose used product as a regular waste, or in accordance with the local waste protocol. The product can also be incinerated, recycled or disposed into a landfill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

Product Symbols



PACKAGING INFO

REF	EAN-13 Barcode	CleanPak™	Carton	Gross Weight
07000	5051400070006	1 pc/bag	25 pcs/carton	~ 4.25 kg

Components/Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
WardPak™	PE	--	√	(√)		√	
Outer Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly recommended (√) = Recommended

ADDITIONAL INFO

- HS Code: 3922900000
- Documentation available upon request. Kindly contact: info@curas.com