

ECLIPSE



MMS
MEDICAL

ECLIPSE



The Eclipse is a rigid manual wheelchair designed for bariatric users. Available in a wide array of sizes, the Eclipse supports up to 1000 lbs. Forward extended casters achieve optimal weight balance on the rear wheels for excellent maneuverability. Overall width is one of the narrowest available allowing for improved door clearance.

A Transportable

Removable wheels, backrest and armrests and front riggings for easy transportation.




B Forward Wheelbase

Allows for more weight to be distributed on the rear wheels, making it easier for foot, hand and attendant propulsion.



C Adjustable Seat Depth and Back Angle

Seat-angle adjustability on chair for comfort and posture control.



D Narrow Overall Width

Compact chair designed to allow easy passage through narrow doorways.



Features

Dynamic back recline option	Removable armrests and front rigging for easy transport
Forward wheelbase for maneuverability	Rigid frame for strength and stiffness
Narrow overall width	Custom seat widths available
Removable back and rear wheels for easy transport	HD Package available for 1000 lbs capacity



SEAT WIDTH
20"-48"



MAX CAPACITY ADJUSTABLE FRAME
600 lbs/272 kg
HD PACKAGE
1000 lbs/453 kg



SEAT DEPTH
16"-24"



BACK HEIGHT
25" & 30"



FRONT SEAT HEIGHT
14.5"-20"



CHAIR WEIGHT
44 lbs/20 kg



BACK ANGLE
85°-100°



OVERALL WIDTH
7.5"