

EDGE 3 *Stretto*

EDGE® 3 STRETTO
shown with iLevel® power seat elevation,
Q-Logic 3 environmental controls, and
Stealth® TRU-Comfort 2 cushion.



The Edge® 3 Stretto features a narrow width for outstanding maneuverability and SRS (Smooth Ride Suspension), making it the perfect indoor and outdoor chair. The Stretto is available with 12.5" or 14" drive wheels and optional 12" of iLevel® power adjustable seat height.



MMS
MEDICAL



SUGAR PLUM

Q
QUANTUM®

Tel: 021 4618000

| www.mmsmedical.ie

Specifications

DRIVE WHEELS	14" or 12.5" pneumatic or flat free tires
CASTER WHEELS	Front: 5" flat free Rear: 6" flat free
MAXIMUM SPEED ¹	Up to 6 mph Up to 3.5 mph w/iLevel activated
GROUND CLEARANCE	2.9"
TURNING RADIUS ²	24.5"
OVERALL LENGTH ²	32.5" without front riggings
BASE WIDTH	21.75" with 14" drive wheels 20.5" with 12.5" drive wheels
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 12-20" Depth: 12-20"
TRU-Balance® 3 Power Tilt	Width: 12-20" Depth: 12-20"
TRU-Balance® 3 Tilt and iLevel®	Width: 12-20" Depth: 12-20"
TRU-Balance® Flex Seating	Width: 15-20" Depth: 15-20" Smaller sizes coming soon
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat ⁵	16-18"
TRU-Balance® 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance® 3 Tilt & iLevel® ⁵	17.5-19.5"
TRU-Balance® Flex Seating	16.5-18.25"
MANUAL RECLINE ⁶	Yes
BATTERY SIZE ³	40 AH (2) or 62 AH (2)
BATTERY WEIGHT ⁷	40 AH: 2 st 2lbs (13.6 kg) Each 62 AH: 2 st.13lbs (19 kg) Each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX 75A Q-Logic NE+ 90A R-Net 70A VR2
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel ⁶
WEIGHT CAPACITY ⁴	21 stone 6 lbs.
BASE WEIGHT	9 stone.

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 21 stone 6 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.
Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability.


Features

- SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- Fender lights and USB charger included
- 10 bold, contemporary colour choices
- Narrow 20.5" base width

Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- Q-Logic 3 EX Drive Controls
- TRU-Balance® 3 Seating
- Trilogly vent tray⁹
- Compact vent tray⁹
- 12" x 12" seat sizes with iLevel (SOF only)

Colours*



* For full range of colours please refer to the Order Form