




Ope-Measure Hoist Scales

Digital Medical Lift Scale gives mobile hoists and ceiling hoist systems the capability of weighing the user.

It enables you to weigh the patients during the moving procedure within just a second. Easy to use, simply press upon the 2 buttons to finish the procedure of weighing and receive the accurate result.

The scale has HOLD and TARE function, and use ordinary AAA batteries. CE-marked and medically approved.



Specifications

Model	Ope-Measure Hoist Scales
Function	Changing Places, Changing Rooms, Assisted Bathrooms
Weight Capacity	300 kg / 400 kg
Scale unit	KG
Division	0.1kg
Outer diameter	30mm
Weight	0.65kg
Power	3 AAA Batteries
Functions	TARE / HOLD
Application	Nursing Homes, Hospitals, Hospices, Domestic Homes



MMS Medical, 51 East Gate Drive, Little Island, Cork T45 EA06
 Tel: 021 4618000 Email: info@mmsmedical.ie www.mmsmedical.ie

Product Specifications Ope-Measure Hoist Scales