



51 East Gate Drive,
Little Island,
Cork, Ireland

T: 021 4618000
E: ppe@mmsmedical.ie



SURGICAL APRON - WATERPROOF & BREATHABLE - WASHABLE AT 95°C

| | | | |
|--|---|--------------------------|----------------------|
| Ref.: 570 020 02 000000 | | TECHNICAL SPECIFICATIONS | |
| COMPOSITION | | UNE EN ISO 1833:2011 | |
| FABRIC | Polyester 100% | BASE: | Polyuretane 100% |
| STRUCTURE | | | |
| FABRIC | Knitted Interlock. Gauge 20. Polyester yarn 75/36/1 Deniers texturized . Total weight 65 g/m ² | | |
| BASE: | Waterproof & breathable membrane 15μ | | |
| TOTAL WEIGHT | 95 g/m ² : 65 g polyester, 19 g polyurethane membrane and 11 g reactive polyurethane glue | | |
| FINISH | Fabric with optical white and Hot Melt lamination with waterproof and breathable membrane | | |
| CHARACTERISTICS | | | |
| DIMENSIONAL STABIL | (Washed at 95°C) Weft-1%. Warp -2% | | UNE EN ISO 5077:2008 |
| WATERPROOFNESS | > 500 cm/H ₂ O | UNE EN ISO 20811:1993 | |
| BREATHABILITY | > 150 g/m ² .24h | BS 7209:1990 | |
| CERTIFICATIONS | | | |
| OEKOT-TEX® Class I ; OEKO-TEX Class II ; 2018OK1127 AITEX / 2002AN4162 AITEX | | | |
| Company certified under ISO 9001:2015 | | | |
| Velamen is member of SEDEX and audited under SMETA protocol ZC1058900 | | | |
| AVAILABLE SIZES: | | Unique size | |
| AVAILABLE COLORS: | | White | |
| INSTRUCCIONS DE MANTENIMENT | | | |
| 30 WASHES AT 95°C GUARANTEED | | | |
| CONFECTION: Overlock elastic safety stitch around the hole for the head | | | |

