

**ROVER<sup>®</sup>**

Made in the Czech republic

**MMS**  
MEDICAL**xRover Reha S****Up to height of 120 cm – up to 44 kg**

CERTIFICATION ČSN EU 1888:2012

[www.mmsmedical.ie](http://www.mmsmedical.ie)

**Variable, ergonomic and comfortable mechanical stroller made from ultra light alloy, especially for children with limited mobility**

Its anatomical shape and conceptual solution ensures for the best possible ergonomic seating position. Pleasant padding satisfies even the most demanding clients with different needs. You can easily adjust the angle of the seat - seating to lying position. The stroller is equipped with anatomic fixation inserts with variable wedge modules, drum brakes with parking lock and many other features for all year-round use in diverse conditions and terrains. This provides the maximum comfort of movement for both the user and the operator.



### Main assets:

- » Support of head and the whole body
- » Safety frame and construction
- » Adjustable handling and seat
- » Extremely low weight
- » Drum brakes with parking lock
- » Comfort and ergonomic seating position
- » Width of stroller – 65 cm
- » The best driving qualities in its class in all terrains

### Optional accessories:

abductwedges, decubitus inserts, additional bags for transport of supporting medical components, protectional elements

### Technical parameters – model S:

Weight (basic equipment): **9,9 kg**  
 Carrying capacity: **44 kg** according to ČSN norm (tested for 70 kg)  
 Adjustable seat: **yes**  
 Length/ width / height: (cm): **118 / 65 / 74**  
 Total height (cm): **108**  
 Total width (cm): **77**  
 Body height: up to **120 cm**

**mms**  
MEDICAL

51 East Gate Drive  
 Little Island, Cork  
 Ireland

T : 021 4618000

E: info@mmsmedical.ie

Volter s.r.o.

Žitná 26

Praha 2

120 00

Czech republic, Europe

Email: info@volter.cz

Tel: +420 737 248 942

**VOLTER**