



A new way to sit!

Zitzi Active is a brand new seat with a special designed back- and seat part which give the User a unique seated position. It is developed to all children and youngsters which have a need for a secure upright seating during activity or while doing work for example at school.

Zitzi Active have a **build-in abductor** which position **the pelvis** to avoid that it rotate backwards and prevent the user from sliding forward.

The forward leaning seat position makes the User have a better **balance and postural control**.

Zitzi Active can with ease be changed from a active seated position to a resting position. In a forward tilted position the User is sitting in a **rider position** (Legs are divided), which urges to stretch up the upper body and makes it harder to collapse. If You tilt the back backwards You create a **resting position without moving the User out from the seat**.

Zitzi Active has a build in movement and a dynamical function which helps the User with reduced muscle strength and / or spasticity.

- Zitzi Active gives the children and youngsters a better sitting where a traditional working chair is not enough.

Zitzi Active facts:

- Available in four standard sizes
- Adjustable backpart
- Several colors available
- Fits to children / youngsters with a length of 100 - 160 cm
- A lot of accessories like knee support, footrest, hand bow, headrest, armrest, belts and so on.
- Fits to different underframes



Underframe Seastar



Underframe Flipper G2

