

# Lallos options


Options are available **in all sizes**

## Comfortable and ergonomic

- Increase the feeling of content
- Helps stabilize the trunk and head in the frontal plane
- Can be used for different children
- Follow the growth

### Front reducer


optional



- Reduces the sitting depth

### Trunk support

optional



## Encompassing and secure


### Activity Table

optional



Standard color by size, customization possible (optional)

- Use as working or play table
- Promotes the erection of the trunk, especially in cases of axial hypotonia
- Increases the feeling of bodily limitations
- Adjustable in height, inclination and has front and side edges
- Material: wood (laminated panels)




### Seat cushion (standard)

- Increases sitting comfort and visual perception with contrast

### Sloping wedge

optional



## Inclusive

- Tilt of 10°
- Increases the comfort of the person by tilting the entire Lalloo backward (the posterior support surface is increased, keeping the trunk and head is easier)
- Standard colors by size
- Anti-slip coating

Suitable for everyone : from baby to adult



# Informations

Laloo is a medical device that offers an alternative to the ordinary sitting position. It's an ergonomic seat in 5 sizes from baby to adult and can be used with legs crossed or outside of the seat. Due to its form, it's encompassing, reassuring and allows the person to be more calm by inviting him to relax in a comfortable seat.

Laloo brings a real inclusive solution and integrates into many worlds such as specialized centers, hospitals, medical offices, waiting rooms, nurseries, schools, relaxation areas, wellness areas, homes...



## Medical Device



Class 1

## Design and Manufacture




French

## Materials



Polyurethane foam (PU)  
PVC cover, Cotton,  
removable


## Warranty



2 years

# Sizes and Colors

Model	Age	Weight	Size	Standards colors
L1	9 - 24 months	7 kg - 11 kg	65 cm - 86 cm	Seat <span style="color: #e91e63;">■</span> - Seat cushion <span style="color: #009688;">■</span>
L2	2 years - 6 years	10 kg - 20 kg	75 cm - 112 cm	Seat <span style="color: #009688;">■</span> - Seat cushion <span style="color: #e91e63;">■</span>
L3	6 years - 12 years	16 kg - 48 kg	100 cm - 158 cm	Seat <span style="color: #ff9800;">■</span> - Seat cushion <span style="color: #fff176;">■</span>
L3.5	12 years - adults	28 kg - 85 kg	150 cm - 180 cm	Seat <span style="color: #fff176;">■</span> - Seat cushion <span style="color: #ff9800;">■</span>
L4	adults	> 80 kg	175 cm - 200 cm	Seat <span style="color: #0070c0;">■</span> - Seat cushion <span style="color: #009688;">■</span>



Color options  
(optional)

