

Permobil

# Panthera Bambino 3



**permobil**

# Empowers the ability to explore

The Panthera Bambino 3 is designed for children who need a robust chair for their active lives to help them gain greater independence. An active, easy-to-handle manual wheelchair for young users aged approximately 4- to 12-years.

## Key benefits:

- Easy to maneuver
- Light weight and easy rolling
- Integrated anti-tip device
- Accessories for added comfort



## Specifications

	Bambino 3	Bambino 3 short
Seat width	24 - 27 - 30 - 33 cm	24 - 27 cm
Seat depth	25 - 30 cm	20 - 25 cm
Seat height	Front 41 cm (22" rear wheels) Back 39 cm (22" rear wheels)	Front 41 cm (22" rear wheels) Back 39 cm (22" rear wheels)
Seat angle	4°	4°
Total length	69 cm	64 cm
Total width	Seat width +28 cm (6° camber)	Seat width +28 cm (6° camber)
Front frame angle	100°	100°
Backrest height	20 - 28 cm or 27 - 35 cm	20 - 28 cm or 27 - 35 cm
Backrest angle	11.5° - (-7.5°)	11.5° - (-7.5°)
Footplate distance	16 - 30 cm	16 - 30 cm
Camber	4° (seat width 27 - 33 cm), 6°, 8°	4° (seat width 27 cm), 6°, 8°
Weight without rear wheels	5.7 - 6.4 kg	5.7 - 6.0 kg
Wheelchair weight	9.8 - 10.4 kg	9.8 - 10 kg
Maximum load	60 kg	60 kg
Material (chassis)	Chrome molybdenum tubing	Chrome molybdenum tubing
Material (upholstery)	Polyurethane covered polyester	Polyurethane covered polyester

Permobil reserves the right to change specifications and design without prior notice. \*More color options available.

Personalise your ride\*

ANODIC  
BLACK

TELE-  
MAGENTA

RUBY  
RED

STING  
YELLOW

MEAN  
GREEN

Actual product colors may vary.