

Permobil

Panthera Micro 3



permobil

Everyday, easy-rolling for children

A light-weight, easy-to-handle manual wheelchair for young users aged approximately 1- to 5-years to use indoors and outdoors on flat surfaces under adult supervision. The Panthera Micro 3 is designed to empower young children to explore the world around them.

Key benefits:

- Easy to push and maneuver
- Robust and lightweight
- Accessories for individual comfort



Specifications

	Micro 3	Micro 3 Long
Seat width	24 - 27 cm	24 - 27 cm
Seat depth	15 - 20 cm	20 - 25 cm
Seat height	Front and rear 40.5 cm (18" rear wheels) Front and rear 33 cm (20" rear wheels)	Front and rear 40.5 cm (18" rear wheels) Front and rear 33 cm (20" rear wheels)
Seat angle	0°	0°
Total length	63 cm	71 cm
Total width	Seat width +22 cm	Seat width +22 cm
Backrest height	22 - 28 cm	20 - 28 cm
Backrest angle	90°	90°
Footplate distance	9,5 - 24,5 cm	9,5 - 24,5cm
Camber	10°	10°
Weight without rear wheels	2.2 - 2.4 kg	2.3 - 2.4 kg
Wheelchair weight	3..4 - 3.5 kg	3.5 - 3.6 kg
Maximum load	30 kg	30 kg
Material (chassis)	Chrome molybdenum tubing	Chrome molybdenum tubing
Material (upholstery)	Polyurethane covered polyester	Polyurethane covered polyester

Personalise your ride*

ANODIC
BLACK

TELE-
MAGENTA

RUBY
RED

STING
YELLOW

MEAN
GREEN

Permobil reserves the right to change specifications and design without prior notice. *More color options available.

Actual product colors may vary.

permobil

permobil.com | info@permobil.co.uk | T. +44 (0)1484 722 888



PANTHERAMICRO3DATASHEETUK221228