

Permobil

# Panthera U3 Light



permobil

# Optimal ergonomics to enable your lifestyle

A low-weight, maneuverable manual wheelchair, the Panthera U3 Light offers stability and comfortable posture support. Designed for optimal ergonomics when seated or driving, this wheelchair can also be easily lifted into a car thanks to its balanced, effortless-grip chassis. Available as U3 Light and U3 Light L.

## Key benefits:

- Ultra-light weight
- Extremely maneuverable
- Easy to transport
- Adaptable to individual needs
- Easy to maintain



## Specifications

	U3 Light	U3 Light Y-front
Seat width	33 - 36 - 39 - 42 - 45 cm	33 - 36 - 39 - 42 - 45 cm
Seat depth	35 - 46 cm	35 - 46 cm
Seat height	Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels)	Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels)
Seat angle	7°	7°
Total length	82 cm	82 cm
Total width	Seat width +21 cm	Seat width +21 cm
Front frame angle	100°	100°
Backrest height	25 - 35 cm	25 - 35 cm
Backrest angle	11.5° - (-4.5°)	11.5° - (-4.5°)
Footplate distance	35 - 43 cm	35 - 43 cm
Camber	2.2°	2.2°
Weight without rear wheels	3.4 - 4.0 kg	3.7 - 3.9 kg
Wheelchair weight	5.7 - 6.3 kg	6.0 - 6.2 kg
Maximum load	100 kg	100 kg
Material (chassis)	Chrome molybdenum tubing	Chrome molybdenum tubing
Material (upholstery)	Polyurethane covered polyester	Polyurethane covered polyester



Personalise your ride\*

ANODIC  
BLACK

PEARL  
WHITE

GREY  
METALLIC

BLUE  
CLASSIC

PINE  
GREEN

Actual product colors may vary.

#### U3 Light Y-front 90°

33 - 36 - 39 - 42 - 45 cm

35 - 46 cm

Front 47 cm (24" rear wheels)  
Back 43 cm (24" rear wheels)

7°

73 cm

Seat width +21 cm

90°

25 - 35 cm

11.5° - (-4.5°)

35 - 43 cm

2.2°

4.0 - 4.2 kg

6.3 - 6.5 kg

100 kg

Chrome molybdenum tubing

Polyurethane covered polyester

#### U3 Light Y-front High 4°

33 - 36 - 39 - 42 - 45 cm

35 - 46 cm

Front 51 cm (24" rear wheels)  
Back 45 cm (24" rear wheels)

4°

82 cm

Seat width +21 cm

100°

25 - 35 cm

11.5° - (-4.5°)

36 - 44 cm

2.2°

3.7 - 3.9 kg

6.0 - 6.2 kg

100 kg

Chrome molybdenum tubing

Polyurethane covered polyester

#### U3 Light Y-front Low

33 - 36 - 39 - 42 - 45 cm

35 - 46cm

Front 44.5 cm (w. 24" rear wheels)  
Back 40.5 cm (w. 24" rear wheels)

7°

82 cm

Seat width +21 cm

100°

25 - 35 cm

11.5° - (-4.5°)

36 - 44 cm

2.2°

3.7 - 3.9 kg

6.0 - 6.2 kg

100 kg

Chrome molybdenum tubing

Polyurethane covered polyester



Enjoy  
the ride

Panthera  
**U3 Light**

For experienced manual wheelchair users who need a low-weight, maneuverable mobility device, the Panthera U3 Light offers stability and comfortable posture support.

Pair it with Mototronik\* power assist device for added convenience.

\*Mototronik - see picture

**permobil**

permobil.com | info@permobil.co.uk | T. +44 (0)1484 722 888



PANTHERAU3LIGHTDATASHEETUK221228