

Universal Deluxe Sling

S4

- Full body supportive sling
- Angled padded leg sections offer full femur & hip support
- Ideal for general transfers
- Small aperture for hygiene purposes
- Colour coded tapes
- Wash temp 85°C / 185°F
- Tumble dry low heat (max 50°C)

Product Codes & Specifications:



	100% Polyester	100% Polyester Mesh	A	B	C	D
XXS	S4AXXSM	S4BXXSM	66 - 76cm	38 - 44cm	51cm	55cm
XS	S4AXSM	S4BXSM	76 - 86cm	44 - 50cm	61cm	60cm
S	S4ASM	S4BSM	86 - 96cm	50 - 57cm	69cm	65cm
M	S4AMED	S4BMED	96 - 108cm	57 - 64cm	76cm	69cm
L	S4ALRG	S4BLRG	108 - 130cm	64 - 69cm	83cm	74cm
XL	S4AXLRG	S4BXLRG	130 - 144cm	69 - 74cm	93cm	79cm
XXL	S4AXXLRG	S4BXXLRG	144 - 160cm	69 - 74cm	103cm	84cm

