

Jazzy Carbon



Budget friendly lightweight powerchair

The Jazzy Carbon Injection is a budget friendly lightweight powerchair weighing only 17.2kg. With it's easy to use controls and folding mechanism it makes the perfect travel companion to meet anyones budget.

Extra lightweight innovative powerchair

Specifications	
Scooter	
Weight Capacity	21 st 4lbs (136kg)
Overall Length	950 mm (37.4")
Overall Width	610 mm (24")
Drive Tyres	Front: 178 mm (7") Solid Rear: 216 mm (8.5") Solid
Ground Clearance	74.9 mm (3 --")
Colour(s)	Black
Mobility	
Turning Radius	900 mm (35.4")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	10 mm
Battery	
Battery Requirements	12 Ah Lithium (flightsafe)
Battery Range	9.3 miles (15 km)
Battery Weight	2.1 kg
Total weight W/O batteries and cushion	17.2 kg


Features




Honeycomb injection technology for extra strength without the cost



Flight safe batteries



Easy-fold mechanism



Easy-to-use controls



REV 01 - 30/06/2025

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®

www.pride-mobility.co.uk