

**Jazzy** Ultra-Lite  
Carbon

Weighing  
only 14.5kg



Colours Available:



# The Ultimate Lightweight Chair

## Travel made easy

This carbon fibre power chair has plenty to offer with its one touch fold release it's never been easier to get moving, whilst weighing only 14.5 kg making it our lightest powered wheelchair yet.

### Specifications

Specifications	
<b>Chair</b>	
Weight Capacity	21 st 5 lb (136 kg)
Opened Dimensions	1000 x 620 x 950 (mm)
Folded Dimensions	620 x 340 x 800 (mm)
Drive Tyres	Front 203 mm (8") Rear 304 mm (12")
Colour(s)	Black
<b>Power</b>	
Turning Radius	900 mm (35.4")
Max Speed	4 mph (6 km/h)
<b>Battery</b>	
Battery Range	12.4 miles (20 km)
Battery Requirements	10 Ah Li (standard) 20 Ah Li (optional)
Battery Weight	10 Ah: 1.4 kg (standard) 20 Ah: 3 kg (optional)
Battery Charger	29.4v 2A
Total Weight Without Batteries	14.5 kg


### Features




Flight safe batteries



Fold away footplate



Compact one lever  
fold function



Stretch mesh backrest for  
lightweight extra comfort



REV 02 - 17/02/2025

**PRIDE**  
MOBILITY PRODUCTS LTD.  
LIVE YOUR BEST®