


# Ramp guides

## System ramps

Fig 1.1

The figure follows the Swedish Institute of Assistive Technology's guidelines for ramp construction.

- The length of the ramp must be at least 12 times the height, however it is recommended that it is 20 times.
- A resting spot is required after 6 metres.



## Portable ramps

Fig 1.2

Recommended models

| Height 10 cm/3,94 in | Height 15 cm/5,90 in | Height 20 cm/7,87 in | Height 30 cm/11,81 in | Height 40 cm/15,75 in | Height 60 cm/23,62 in |
|----------------------|----------------------|----------------------|-----------------------|-----------------------|-----------------------|
| X050                 | IRC075               | X100                 | F150                  | T200                  | T300                  |
| IRP050               |                      | F 110                | IRP120                | T 200-3               | FTW300                |
|                      |                      | XW100                | IRP150                | FW200                 | FT300                 |
|                      |                      | IRP090               | IRP170                | FTW200                | VR28                  |
|                      |                      | IRP120               | T120                  | IRP170                |                       |
|                      |                      | T090                 | IRC120                | IRP200                |                       |
|                      |                      | VR10                 | IRC150                | FT 200                |                       |
|                      |                      | VR10-1               | VR16                  | IRC200                |                       |
|                      |                      |                      |                       | VR20                  |                       |
|                      |                      |                      |                       | VR25                  |                       |