



**SUPREME PERFORMANCE
AND EXHILARATING STYLE**



Distributed by:

MMS
MEDICAL



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

Pride®
Mobility Products Ltd.
— LIVE YOUR BEST® —



CORK Tel: 021 461 8000 | Fax: 021 461 8099

DUBLIN Tel: 01 450 0582 Fax: 01 456 5889

Email: info@mmsmedical.ie

Web: www.mmsmedical.ie

(Version 1: May 2017)

Features and Specification

Introduction

The sleek, sporty 6 mph Colt Deluxe mobility scooter from the 'Pride® Scooter' range delivers a smooth performance. With all-new advanced standard features including front and rear suspension, high-visibility automatic brake lights, feather-touch disassembly, wraparound delta tiller and LED kerb light.



Available on



Available Colours

The Colt Deluxe comes in two optional colours:














Red



Titanium Grey



Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
 28 Stone 8lbs. (181 kg)	 6.25 Mph (10 km/h)	 5.7 cm (2.25 inch's)	 133.35 cm (52.5 inch's)	 119.4 cm (47 inch's)	 56 cm (22 inch's)	 Front: 25.4cm Rear: 25.4 cm	 10.7 Stone 68 k.g.	 29.48 k.g. Front Section	 Hi-Back Seat 19 k.g	 x2 40-45 Ah Battery

Warranty

1. Life time limited warranty on frame
2. Two year limited warranty on drive train, motor and brake
3. One-year limited warranty Electronic Components, battery and Charger.

Notes

Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. Battery weight may vary based on manufacturer.



MMS Medical
51 East Gate Drive,
Little Island,
Cork.

MMS Medical
Unit 81
Western Parkway Business Park,
Lower Ballymount Road,
Dublin 12

CORK Tel: 021 461 8000 | Fax: 021 461 8099
DUBLIN Tel: 01 450 0582 Fax: 01 456 5889
Email: info@mmsmedical.ie
Web: www.mmsmedical.ie

