



Care and treatment table MONA
The optimum, mobile solution for demanding care

Care and treatment table MONA

the ideal solution for demanding care ...





Continuous, vertical height adjustment for simple transfer and optimum working height ...

A large, high-quality support pad forms the basis for safe and convenient treatment of people in need of assistance. The rounded ends of the table are padded and prevent injuries, even to very susceptible skin. Even if the person is very restless, the sturdy scissors height adjustment ensures that neither the support pad nor the treatment table can wobble or drop.

The treatment table can be furnished with a folding side rail on each side, which provides additional safety. A sturdy metal fastener makes sure that the side rails cannot fold down, even under strong pressure. Continuous height adjustment (optionally electrical or hydraulic) allows the optimum work and transfer height to be set. The base frame can be accessed underneath from all sides and is equipped with four wheels with steering and brakes.



Technical data

Pad length (optional) - special sizes on request	1.500 mm to 2.000 mm
Pad width (optional) - special sizes on request	650 mm / 800 mm / 1.000 mm
Lifting height (hydraulic height adjustment)	490 - 955 mm
Lifting height (electrical height adjustment)	465 - 1.005 mm
Max. load	150 kg (Version EL), 135 kg (Version Hyd)

Accessories / options

Roll-off side rails
Protective pad for side rails
Head section with notch adjustment



MMS Medical Ltd.
51 East Gate Drive
Little Island,
Cork.

Tel: 021 4618000
Fax: 021 4618099
E: info@mmsmedical.ie

MMS Medical Ltd.
Unit 81, Western Parkway Business Park,
Lower Ballymount Road,
Dublin 12.

Tel: 01 4500582
Fax: 01 4565889
www.mmsmedical.ie