

PRODUCT DATASHEET

Patient Wipes, 28x30cm Airlaid, 50gsm, Maceratable, SOFT

Brand	: CURAS®
Category	: Nursing
Range	: Wipes
Product Group	: Dry Wipes
REF	: 02090



PRODUCT DESCRIPTION

Disposable patient wipes made of highly absorbent, maceratable pulp material.

INTENDED USE

Patient wipes are suitable for professional washing of incontinent, bedridden or dependent individuals in hospitals, nursing homes and home healthcare. They are also helpful for outpatient care and maternity. Use with water, soap or lotion, as required. The highly absorbent wipes also help to dry patients quickly and easily. Single use only. Patient wipes can be macerated after use.

FEATURES & BENEFITS

- Cost-efficient, hygienic and an effective infections control measure.
- Colour-free and fragrance-free material.
- Made from formaldehyde-free and totally chlorine-free (TCF) material.
- Versatile, can be used in conjunction with any soap, lotion, or water.
- Highly absorbent, suitable for drying off patients.
- Disposable for reduced risk of cross-contamination and Hospital-Acquired Infections.
- Product is low lint, which means there is less residual fibre after use.
- Macerator friendly, can be easily disposed by maceration.
- Each piece is Z+1/2 folded for easy dispensing
- Inner bag protects product from dust, contamination and liquids at all times.
- Centre perforation on bag allows each wipe to be dispensed easily without full exposure.

PRODUCT SPECIFICATIONS

REF	Gram per Square Meter	Dimensions (L x W)	Average Weight/Pc	Packing
02090	50 g/m ² ± 5 g/m ²	28 x 30 cm (~11 x 12 in)	4.2 g ± 0.4 g	Inner Bag

Composition

- 100% totally chlorine-free (TCF) fluff pulp.
- Colour-free and fragrance-free.
- Formaldehyde-free.

Country of Origin

- China

Sterilization

- Non-sterile

Shelf Life

- 8 years

Selected Standards & Guidelines

- ISO 9001:2008
- GB/T 2912.1-2009 (formaldehyde detection in textiles)

Recommended Storage Condition

- Store at room temperature (5°C-30°C/41°F-170.6°F) and avoid exposure to elevated temperatures.
- Keep in a cool and dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- Dispose of used product into a macerator, or in accordance with local/clinical waste protocol. If a macerator is not available, incinerate or land-fill.
- Do not dispose of used product into the toilet bowl.

Product Symbols



curas.com/Symbols

PACKAGING INFO

REF	EAN-13 Barcode	Inner Bag	Carton	Gross Weight
02090	5051400020902	75 pcs / bag	2,250 pcs / carton (30 inner bags)	~11 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
Inner Bag	PE	--	√	(√)		√	
Outer Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly Recommended (√) = Recommended

ADDITIONAL INFO

- HS Code: 4015110000
- Documentation available upon request. Kindly contact: info@curas.com