

PRODUCT DATASHEET

Disposable Mid-Stream Urine Sampler, 1L, WardPak™

Brand	: CURAS®
Category	: Nursing
Range	: Hygiene Pulp™
Product Group	: Human Waste Collection
REF	: 08040



PRODUCT DESCRIPTION

Disposable Mid-Stream Urine Sampler made from hygienic pulp material with water resistance capability. Non-sterile. Suitable for disposing in macerator.

INTENDED USE

Disposable Mid-Stream Urine Sampler is used for collecting fresh urine sample for the purpose of testing. Suitable for single use only.

FEATURES & BENEFITS

- Disposable and single use can reduce risk of cross-contamination and Hospital-Acquired Infections.
- Flat base design of the bowl allows stable placement.
- Product with thick base to maintain material strength during usage and transportation.
- Available with wide edge for easy holding and handling when content is full.
- Product is stackable for easy storage and reduces storage space.
- Graduation marks provide the convenience of measuring the content, even to small amount, i.e. 100ml.
- V-shaped spout ensures easy & accurate transfer of the content to other containers, without spillage.
- Shape that fit perfectly on commode chair and toilet bowls.

PRODUCT SPECIFICATIONS

REF	Capacity	Bowl Dimensions (L x W x H)	Average Weight / Pc	Inner bag
08040	1L	32.63(L) x 20.6(W) x 9.75(H) cm	38 – 40 g	Wardpak™

Composition

- Made of 100% recycled newsprint paper.
- Water-based fluid resistant compound.

Country of Origin

- China

Sterilization

- Non-Sterile.

Shelf Life

- 8 years

Selected Standards & Guidelines

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- PAS29:1999 Fluid resistance test for a minimum period of 4 hours (water at 35°C ±3°C).
- PAS29:1999 Maceration test with a hospital sluice room type macerator, normal cycle.

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- Dispose of used product into a pulp macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.

Product Symbols

curas.com/Infection-Prevention



curas.com/Symbols

PACKAGING INFO

REF	EAN-13 Barcode	Wardpak™	Carton	Gross Weight
08040	5051400080401	25 pcs / bag	250 pcs / carton	~ 11 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
Inner Bag	LDPE bag	--	√	(√)		√	
Outer Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly recommended (√) = Recommended

ADDITIONAL INFO

- Curas® Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com