

PRODUCT DATASHEET

Holder for Urine Bottle, Universal, Stainless Steel, Autoclavable, Ø40mm

| | |
|---------------|---------------------|
| Brand | : CURAS® |
| Category | : Nursing |
| Range | : Hangers & Holders |
| Product Group | : Hangers & Holders |
| REF | : 09411 |



PRODUCT DESCRIPTION

Holder for urine bottle, made of stainless steel rod. Autoclavable.

INTENDED USE

The stainless steel holder for urine bottle is designed to hold and fit any standard male urine bottles, including Curas pulp urine bottles. Easy to use and autoclavable.

FEATURES & BENEFITS

- Made of genuine stainless steel 304, that is strong and resistant to corrosion.
- Light weighted holder with smooth surface that allows easy cleaning, disinfecting and autoclaving. Thus, conform to strict hygiene requirement.
- Suitable for holding any standard pulp and plastic male urinals, including Curas pulp urine bottles (REF 08050/56/57/58/59/60).
- Convenient and flexible for attaching to any standard bed frame.
- Designed to hold urinal in an upright position without risk of slippage.
- Designed with rounded edge to ensure user's safety and for better holding of the product.
- Packed individually in inner bag for easy dispensing and protection against dust.
- Stackable, and thus reduces storage space.

PRODUCT SPECIFICATIONS

| REF | Dimensions (L x W x H) | Average Weight/Pc |
|-------|--------------------------------|-------------------|
| 09411 | 155 (L) x 135 (W) x 200 (H) mm | ~147g |

Composition

- Made of Stainless Steel 304
- Finishing by Electropolishing
- Non-magnetic

Country of Origin

- China

Sterilization

- Non-sterile
- Can be sterilized by autoclave.

Shelf Life

- 15 years

Selected Standards & Guidelines

- N/A

Recommended Storage Condition

- Store at room temperature (10°C-25°C / 50°F-77°F).
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- Disposal of worn out product to applicable point for recycling of stainless steel product.

Product Symbols

curas.com/Infection-Prevention



curas.com/Symbols

PACKAGING INFO

| REF | EAN-13 Barcode | Inner Bag | Inner Box | Carton | Gross Weight |
|-------|----------------|---------------|-------------|-----------------------|--------------|
| 09411 | 5051400094118 | 1 hanger/ bag | 12 pcs/ box | 36 pcs/ carton (3x12) | ~7kg |

| Components / Material Properties | | | Recommendations on Waste Handling | | | |
|----------------------------------|-------------------------|--------------|-----------------------------------|----------|---------|-------|
| Packaging | Material | Contains PVC | Incinerate | Landfill | Recycle | Reuse |
| Inner Bag | Transparent bag | -- | √ | (√) | √ | |
| Inner Box | White Corrugated Box | -- | (√) | (√) | √ | |
| Carton | Brown Corrugated Carton | -- | (√) | (√) | √ | |
| Transport Pallet | Wood | -- | | | (√) | √ |

√ = Highly recommended

(√) = Recommended

ADDITIONAL INFO

- HS Code: 7326909890
- Documentation available upon request. Kindly contact: info@curas.com