

Q6 edge[®] 2.0



A revolutionary redesign to the industry renowned Q6 Edge[®].

The new Q6 Edge[®] 2.0 features versatile powered solutions, which can be configured with the highly adaptable TRU-Balance[®] 3 Seating System. Meeting the needs of the most active user.

UNIQUE TECHNOLOGY

to the industry.

iLevel™ | TRU-Balance[®] Seating System
Q-Logic 2 | TB3 Flip-up armrests

Distributed by:

MMS
MEDICAL



QUANTUM[®]
The Rehab Power Chair Solutions Company



CRASH
TESTED



@QUANTUMREHABUK

CORK Tel: 021 461 8000 | Fax: 021 461 8099 | DUBLIN Tel: 01 450 0582 Fax: 01 456 5889
Email: info@mmsmedical.ie | Web: www.mmsmedical.ie



ATX® adjustable suspension

Provides increased ride quality, obstacle climbing and stability
Semi-independent articulating rear castor beam.

55 Ah Batteries

More miles per charge! 55 Ah batteries allow for increased range, for the most active user. Easy access through front shroud.

Improved Aesthetics

Updated front & rear shroud design.
New "snap on" plastic chassis shrouds

Turning Radius

Access smaller spaces easier with a tight turning radius
of 20.5" (52.07 cm)

Motor Package

2-Pole advanced responsive motors
4 Mph or 6 Mph motor options available

Measurements

14" Drive wheel, 6" low impact solid castor wheels front and rear.
Length: 35.5" (Without legrests) | Width: 24.4" |
Ground Clearance: 2.7"

Seating Solutions



Synergy Seating

- 10" x 10" to 22" x 22" seat size options
- Visco padded or tension adjustable backrest options
- Durable and easily adjustable
- Easy centre of gravity adjustment
- Accepts a wide variety of after market seats and backs



TRU-Balance® 3 Seating

- 12" x 12" to 22" x 22" seat size options
- Make fine adjustments while client is in the chair
- Accepts a wide variety of after market seats and backs
- Endomorph and mesomorph seating, allows for independent adjustment of seat width and backrest width
- Premium accessories available

Power Options



TRU-Balance® 3 Power Tilt

- 0-50° of tilt (Anterior tilt available)
- Lower standard seat-to-floor height
- Easy adjustment and durability
- Maintains centre of gravity shift



TRU-Balance® 3 Power Recline

- 85° - 175° of recline
- Adjustable shear reduction
- Compatible with after-market backrest interface

Powered Elevating leg rests

- Articulating foot platform
- Articulating swing-away leg rests
- Can be operated in conjunction with TB3 recline




@sophieslifeonilevel

10" Power Adjustable Seat Lift

- Lower standard seat-to-floor height
- iLevel™ allows 10" of seat elevation while at walking speed (3.5 mph.)
- iLevel is retrofittable on the Q6 Edge 2.0



Additional Features



TRU-Balance® 3 Cantilever

Flip-up armrests

- Quick height adjustment
- Can be internally or externally rotated for optimal positioning
- Multiple armrest pad options



TRU-Balance® 3 Centre

mount foot platform

- Fixed or powered articulating
- Modular footplate extension kits
- Multi-adjustable calf pads
- Paediatric & adult footplate options



Swing-away legrests

- 60° - 90° standard swing-away
- Paediatric & adult footplate options
- Powered Articulating options
- Manual elevating



Attendant/Dual control options

- Full proportional joystick
- Quick and easy to use control interface
- Control the chair at the attendants speed




Body positioning

- A full range of multi adjustable swing-away laterals, offering advanced positioning of the trunk
- Designed to easily allow head to toe positioning using any of the Stealth positioning accessories



After-market backrest interface

- Allows for an unprecedented level of compatibility between TRU-Balance 3 back rest and most after-market backrests
- Retains full functionality, including shear reduction



Premium accessories

- Transfer bars
- Mobile phone holder
- USB Port (available with Q-Logic 2)
- Clothes guard
- Additional accessories available



Full range of control devices available

- Switches
- Sensors
- Proportionals
- i-Drive™



OPT: Extended Traction Tyres⁴



Quantum Q6 Edge 2.0 | Colour Options & Tech Spec

Colour Options



Red Corvette



Back in Black



Pretty in Pink



Lime Light



Silver Bullet



Orange Crush



Electric Blue



White Satin



Deep Purple



Grey Street



Midnight Blue



Yellow Blaze

Technical Specifications

Drive Wheels (cm) 35.56 cm pneumatic or solid	Caster Wheels (cm) Front and Rear: 15.24 cm	Max Speed (mph) 4mph or 6 mph with lights ¹	Ground Clearance (cm) Battery Tray: 6.9 cm	Turning Radius (cm) 52.1 cm	Overall Length (cm) 90.2 cm	Overall Width (cm) 62.2 cm	Seating Options (inches) Synergy and TB3 Seating	Batteries x 2 50 - 55 AH Deep Cycle	Available Electronics VR2/RNet NE+/QLogic 2	Weight Capacity (stone) 21 Stone 6 lbs (137k.g)	
Seat to ground height (") 16.5"-19.5"	Base Weight (kg) 63 kg	Battery Weight Each (kg) 17.24 kg	Warranty <ol style="list-style-type: none"> Lifetime limited warranty on frame. 1-year limited warranty on electronics. 18-month warranty on drive control system. 2-year limited warranty on drivetrain. 				Notes <ol style="list-style-type: none"> Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. Battery weight may vary based on manufacturer. Not available with iLevel Seat lift 				

Get Social with Quantum



FOLLOW THE LIVES OF OUR QUANTUM STARS, AND GET THEIR PERSPECTIVE OF THE Q6 EDGE 2.0

www.quantumrehab.co.uk/quantum-stars



Ensure you have the most up to date information about Quantum Rehab by subscribing to our social media pages

Distributed by:



MMS Medical
51 East Gate Drive,
Little Island,
Cork.

MMS Medical
Unit 81
Western Parkway Business Park,
Lower Ballymount Road,
Dublin 12

Further information can be obtained through our web site, or prescription forms. All information is subject to change without notification. Please consult Quantum with any queries you may have.