



1 Crash tested - and approved!

Seabass PRO is crash tested and approved according to ISO 7176-19 and ISO 10542.




SWEDISH QUALITY AND DESIGN



2 Seat depth adjustment

Finally it's possible to adjust the seat depth both forwards and backwards so that you can easily find the perfect balance between this chassi and any number of seating systems and/or custom made seats.



3 Tilt in space

Unique for Seabass PRO is its "tilt in space" feature, allowing the seat to be tilted without compromising the balance and stability between the seat and chassi. This way the center of gravity remains unchanged so that the seat system is always perfectly balanced regardless of the degree of tilt.



4 Wheels

In order to fulfill our customers' wishes, rear wheels are available in 5 different sizes. Choose between 12", 16", 20", 22", and 24" wheels. All of the rear wheels are equipped with quick-release function and are available with either air-filled or solid tires.

Picture is showing
 Seabass Pro size 2 with 22" quick-release wheels and multi-adjustable footrest



Driving bow
 Three different models and sizes

Driving bow
 You can easily adjust both height and angle

Wheels
 Chose between 5 different sizes

Strong
 Firm construction

Quick-release
 For easy mounting

Seabass Pro
 Strong, balanced, stable and cool design

 **SWEDISH QUALITY AND DESIGN**

Seabass PRO Quick facts

Seabass Pro with its center of gravity on the right place. Versatile wheel location for a huge variety and wheel sizes. Whatever seating solution chosen this guarantees perfect balance for every fitting.

Sizes	3 st	Seat height	45 cm
Wheels	12-24"	Color	Black
User weight	135 kg	Weight	19 kg
Tilt	30°	Crash-tested	Approved